

# DERMAFLEX DFS2

DERMAFLEX DFS2 is an anionic neutral syntan based on a dicyandiamide resin.

## TYPICAL ANALYSIS

Appearance	Fine white powder
Solid content	> 92 %
pH (analytical solution)	$8.0 \pm 1.0$
Lightfastness:	Excellent
Heatfastness:	Excellent

## PROPERTIES

DERMAFLEX DFS2 is recommended as a retannage for neutralised chrome stock. It is fixed under acidic conditions and has a selective filling action which improves the more open structured areas of leather.

This results in a leather with a full feel, tight flexible grain and uniform softness. It helps to promote level dyeing and even uptake of fatliquors.

DERMAFLEX DFS2 may be used on bovine or small skin wet blue for full grain leathers. Moreover it will assist in even buffing for corrected grain or nubuck leathers, in conjunction with other retanning agents.

DERMAFLEX DFS2 presents no problem under normal storage conditions and has an indefinite shelf life. Some caking may occur if the product is stored in warm damp conditions but this does not detract from its efficiency.

## APPLICATION

Recommended offer from 2 to 5 %. For corrected grain or nubuck leathers an offer of 3 to 5 % is suggested, in combination with other retanning agents.

When correctly stored, it is stable 36 months.

This sheet is a synthetic illustration of the product to enable the user to employ it at the best.

Properties and values are not, therefore, sales specifications.

Taking into account the several possibilities of application of this product, no warranty of the results is possible.

**Expiry date refers to the product well closed in the original container, in fresh and dry area.** Some particularly sensitive products, when open, should be used rapidly.

**Stir well before using the product, close tightly the container to preserve the properties.**

Our technical service is at your disposal for any recommendation.